



Cylinder Storage & Transportation

Cylinder Trolleys 88

Cylinder Wall Brackets, Holders & Plates 91

Cylinder Trolley F/G Size

Features

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

Overall Dimensions (w x d x h)

- 510 x 405 x 1160mm

Standard Colours

- Paint - Grey



CYT/FG1

Item	Description	Weight
CYT/FG1	F/G Size	10.0kg

Cylinder Trolley with Stabilisers F/G Size

Features

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

Overall Dimensions (w x d x h)

- 520 x 700 x 1160mm

Standard Colours

- Paint - Grey



CYT/FG1/ST

Item	Description	Weight
CYT/FG1/ST	F/G Size	14.0kg

Cylinder Trolley with Stabilisers F/G Size Double

Features

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinders
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 200mm anti-static wheels, 100mm swivel castors

Overall Dimensions (w x d x h)

- 780 x 740 x 1145mm

Standard Colours

- Paint - Grey



CYT/FG2/ST

Item	Description	Weight
CYT/FG2/ST	F/G Size Double	20.5kg

Cylinder Trolley J Size

Features

- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter 245mm, height 1550mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

Overall Dimensions (w x d x h)

- 510 x 475 x 1280mm

Standard Colours

- Paint - Grey



CYT/J1

Item	Description	Weight
CYT/J1	J Size	10.5kg

Cylinder Trolley with Stabilisers J Size

Features

- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter, max diameter 245mm, height 1550mm
- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

Overall Dimensions (w x d x h) Standard Colours

- 510 x 740 x 1275mm

- Paint - Grey



CYT/J1/ST

Item	Description	Weight
CYT/J1/ST	J Size	14.8kg

Bulk Transportation Trolley D & E Cylinders, 12 x Capacity

Features

- Compatible with D, AD, CD, DD, PD, RD, ZD & E (24 cu ft) cylinders (12 x capacity) max diameter 105mm, height 870mm
- One piece construction
- Cylinders locate into individual sections, supported on a non slip rubber base
- 125mm swivel castors, 2 off braking

Overall Dimensions (w x d x h) Standard Colours

- 460 x 590 x 990mm

- Paint - Grey



CYT/DE12

Item	Description	Weight
CYT/DE12	D & E Cylinders, 12 x Capacity	23kg

Cylinder Trolley E Size

Features

- Compatible with E (24 cu ft) cylinders, max diameter 115mm, height 900mm
- One piece construction
- Recessed base plate and drop over chain to secure the cylinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

Overall Dimensions (w x d x h)

- 320 x 235 x 980mm

Standard Colours

- Paint - Grey



CYT/E1

Item	Description	Weight
CYT/E1	E Size	4.3kg

Cylinder Trolley CD/D Size

Features

- Compatible with CD & D cylinders, max diameter 115mm, height 540mm
- One piece construction
- Recessed base plate and drop over chain to secure the cylinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

Overall Dimensions (w x d x h)

- 320 x 235 x 980mm

Standard Colours

- Paint - Grey



CYT/CD1

Item	Description	Weight
CYT/CD1	CD/D Size	4.3kg

Cylinder Trolley AF Size

Features

- Compatible with AF, DF & ZH cylinders, max diameter 178mm, height 900mm
- One piece construction
- Recessed base plate and drop over chain to secure the cylinder
- 100mm anti-static wheels

Overall Dimensions (w x d x h)

- 390 x 310 x 980mm

Standard Colours

- Paint - Grey



CYT/AF1

Item	Description	Weight
CYT/AF1	AF Size	4.9kg


Wall Brackets Cylinder

Features

- Max cylinder diameters
 - D/E - 112mm, F - 148mm, G - 185mm, J - 235mm
- One piece construction
- Drop over chains to secure the cylinder
- Holes for wall mounting (fixings not provided)
- Should be positioned on the wall at approx half cylinder height

Standard Colours

- Paint - Grey



Item	Description	Dimensions (w x d x h)	Weight
CYB/DE2	2 x D/E Cylinders	359 x 107 x 88mm	1.1kg
CYB/F2	2 x F Cylinders	435 x 126 x 88mm	1.4kg
CYB/G2	2 x G Cylinders	511 x 145 x 88mm	1.7kg
CYB/J2	2 x J Cylinders	613 x 170 x 88mm	2.1kg

Cylinder Holder D & E Size, Vertical

Features

- Compatible with...
 - D Size - D, AD, CD, DD, PD, RD & ZD cylinders (max 112mm diameter)
 - E Size - E (max 112mm diameter)
- One piece construction
- Carry handle and extended hook
- Suitable for beds with 25mm wide side/end members
- Holder can be hung and/or bolted to 25mm section bed-frame

Overall Dimensions (w x d x h)

- CYH/D1 - 200 x 160 x 290mm
- CYH/E1 - 200 x 160 x 595mm

Standard Colours

- Paint - Grey



CYH/D1

CYH/E1

Item	Description	Weight
CYH/D1	D Size, Vertical	1.4kg
CYH/E1	E Size, Vertical	2.5kg

Cylinder Holder Wall Plate

Features

- Painted stainless steel wall plate to support cylinder holders
- Compatible with CYH/D1, CYH/E1, CYHP/D1 & CYHP/E1

Overall Dimensions (w x d x h)

- 240 x 47 x 254mm

Standard Colours

- Paint - Grey



CYWP/1

CYWP/1 Shown with CYH/D1

Item	Description	Weight
CYWP/1	Wall Plate	2.1kg