



Dressing / Instrument Trolleys

Removable Shelves 102

Fixed Shelves 104

Drawer Units 107

Options 108



Removable Shelves, Stainless Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or a shelf with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Options

- C001 - 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 - 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 - 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 - 75mm, Pressed Steel Castors (4 x braking)
- C014 - 75mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- Drawer Units - Range of sizes / configurations (see page 107)
- DP/C - Infusion Pole & Bracket
- E/SOCKET/4 - Four Gang Socket (requires an additional drawer)
- GR - Guard Rail
- ZP005 - Ring Buffers

Size	Flat Shelves	Reversible Shelves	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSR/450	DTSR/450/SR	450 x 450mm	470 x 470 x 890mm	11.5kg
Medium	DTSR/600	DTSR/600/SR	600 x 450mm	620 x 470 x 890mm	16.0kg
Large	DTSR/750	DTSR/750/SR	750 x 450mm	770 x 470 x 890mm	18.0kg
X-Large	DTSR/900	DTSR/900/SR	900 x 450mm	920 x 470 x 890mm	20.0kg

Removable Shelves, Mild Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded mild steel frame
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or a shelf with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Standard Colours

- Paint - White

Options

- C001 - 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 - 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 - 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 - 75mm, Pressed Steel Castors (4 x braking)
- C014 - 75mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- Drawer Units - Range of sizes / configurations (see page 107)
- DP/C - Infusion Pole & Bracket
- E/SOCKET/4 - Four Gang Socket (requires an additional drawer)
- ZP005 - Ring Buffers



Size	Flat Shelves	Reversible Shelves	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTMR/450	DTMR/450/SR	450 x 450mm	470 x 470 x 890mm	11.5kg
Medium	DTMR/600	DTMR/600/SR	600 x 450mm	620 x 470 x 890mm	16.0kg



Fixed Shelves, Sides Up, Stainless Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Options

- C001 - 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 - 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 - 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 - 75mm, Pressed Steel Castors (4 x braking)
- C014 - 75mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- Drawer Units - Range of sizes / configurations (see page 107)
- DP/C - Infusion Pole & Bracket
- E/SOCKET/4 - Four Gang Socket (requires an additional drawer)
- ZP005 - Ring Buffers

Size	Sides Up	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/450	450 x 450mm	470 x 470 x 890mm	11.5kg
Medium	DTSF/600	600 x 450mm	620 x 470 x 890mm	16.0kg
Large	DTSF/750	750 x 450mm	770 x 470 x 890mm	18.0kg
X-Large	DTSF/900	900 x 450mm	920 x 470 x 890mm	20.0kg

Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Available with one, two or three shelves
 - One shelf - 870mm (lower central tie tube minimises surface area for cleaning, ideal for use in theatres)
 - Two shelf - 405 & 870mm
 - Three shelf - 255, 560 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Options

- C001 - 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 - 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 - 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 - 75mm, Pressed Steel Castors (4 x braking)
- C014 - 75mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- Drawer Units - Range of sizes / configurations (see page 107)
- DP/C - Infusion Pole & Bracket
- E/SOCKET/4 - Four Gang Socket (requires an additional drawer)
- ZP005 - Ring Buffers

One Shelf

Size	Sides Down	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/450/1/SD	450 x 450mm	470 x 470 x 870mm	9.9kg
Medium	DTSF/600/1/SD	600 x 450mm	620 x 470 x 870mm	11.7kg
Large	DTSF/750/1/SD	750 x 450mm	770 x 470 x 870mm	14.5kg
X-Large	DTSF/900/1/SD	900 x 450mm	920 x 470 x 870mm	16.8kg

Two Shelf

Size	Sides Down	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/450/2/SD	450 x 450mm	470 x 470 x 870mm	12.0kg
Medium	DTSF/600/2/SD	600 x 450mm	620 x 470 x 870mm	14.5kg
Large	DTSF/750/2/SD	750 x 450mm	770 x 470 x 870mm	18.0kg
X-Large	DTSF/900/2/SD	900 x 450mm	920 x 470 x 870mm	21.0kg

Three Shelf

Size	Sides Down	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/450/3/SD	450 x 450mm	470 x 470 x 870mm	15.3kg
Medium	DTSF/600/3/SD	600 x 450mm	620 x 470 x 870mm	19.9kg
Large	DTSF/750/3/SD	750 x 450mm	770 x 470 x 870mm	23.5kg
X-Large	DTSF/900/3/SD	900 x 450mm	920 x 470 x 870mm	27.6kg





Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Available with one, two or three shelves
- One shelf - 870mm (lower central tie tube minimises surface area for cleaning, ideal for use in theatres)
- Two shelf - 405 & 870mm
- Three shelf - 255, 560 & 870mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel anti-static non-marking castors

Options

- C002 - 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 - 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C015 - 100mm, Pressed Steel Castors (4 x braking)
- DP/C - Infusion Pole & Bracket
- ZP005 - Ring Buffers

One Shelf

Size	Sides Down	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/66/1/SD	600 x 600mm	620 x 620 x 870mm	14.0kg
Medium	DTSF/96/1/SD	900 x 600mm	920 x 620 x 870mm	18.0kg
Large	DTSF/126/1/SD	1200 x 600mm	1220 x 620 x 870mm	20.4kg

Two Shelf

Size	Sides Down	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/66/2/SD	600 x 600mm	620 x 620 x 870mm	17.5kg
Medium	DTSF/96/2/SD	900 x 600mm	920 x 620 x 870mm	24.0kg
Large	DTSF/126/2/SD	1200 x 600mm	1220 x 620 x 870mm	29.6kg

Three Shelf

Size	Sides Down	Shelf Size (w x d)	Overall Dimensions (w x d x h)	Weight
Small	DTSF/66/3/SD	600 x 600mm	620 x 620 x 870mm	22.8kg
Medium	DTSF/96/3/SD	900 x 600mm	920 x 620 x 870mm	32.7kg
Large	DTSF/126/3/SD	1200 x 600mm	1220 x 620 x 870mm	41.2kg

Drawer Units

Features

- Compatible with removable & fixed shelf ranges, with exception of 600mm deep trolleys
- Available in three widths...
 - DRW260 - Set of two drawers
 - DRW350 - Single drawer
 - DRW650 - Single drawer
- Maximum capacity 10.0kg
- Factory fitted & must be ordered at the same time as the trolley

Internal Dimensions (w x d x h)

- DRW260 - 230 x 410 x 100mm
- DRW350 - 350 x 410 x 100mm
- DRW650 - 650 x 410 x 100mm

Options

- E/SOCKET/4 - Four Gang Socket



DTSR/750 with DRW260R/DRW260

Removable Shelf Trolleys

Trolley Width	One Tier Combination Code			Two Tier Combination Code		
	DRW260/R	DRW350/R	DRW650/R	DRW260/R DRW260	DRW350/R DRW350	DRW650/R DRW650
450mm						
600mm	✓	✓		✓	✓	
750mm	✓	✓	✓	✓	✓	✓
900mm	✓	✓	✓	✓	✓	✓

Fixed Shelf Trolleys

Trolley Width	DRW260/F	DRW350/F	DRW650/F	DRW260/F DRW260	DRW350/F DRW350	DRW650/F DRW650
	450mm					
600mm	✓	✓		✓	✓	
750mm	✓	✓	✓	✓	✓	✓
900mm	✓	✓	✓	✓	✓	✓

Castors

Features

- Polymer & pressed steel castors
- Anti-static & non-marking
- Supplied as a set of 4 castors
- 80mm polymer castors increase the height of the trolley by approx. 5mm
- 100mm polymer castors increase the height of the trolley by approx. 30mm
- 100mm steel castors increase the height of the trolley by approx. 30mm
- Factory fitted



Item	Description
C001	Castor - 80mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4)
C002	Castor - 100mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4)
C003	Castor - 100mm, Pressed Steel, 2 x Braking / 2 x Non-Braking (set of 4)
C013	Castor - 75mm, Pressed Steel (4 x braking)
C014	Castor - 75mm, Pressed Steel (2 x braking / 2 x non-braking)
C015	Castor - 100mm, Pressed Steel (4 x braking)

Four Gang Socket

Features

- BS EN 60950:2006 Information technology equipment. Safety. General requirements
- Switched & fused four gang socket
- Supplied with 3 metres of beige cable & fitted with a UK 3 pin plug
- Only compatible when fitted in conjunction with an optional drawer
- Fitted to the rear of the drawer, complete with a cleat to secure the cable when not in use
- Factory fitted



Item	Description	Weight
E/SOCKET/4	Four Gang Socket	1.4kg

Ring Buffers

Features

- Set of four ring buffers
- Offers protection to both the trolley & fabric of the building
- Increases the length & width of trolley by approx. 80mm
- Factory fitted



Item	Description	Weight
ZP005	Ring Buffers (set of four)	0.3kg

Reversible Shelves

Features

- Allows the end user to choose between a flat shelf or one with three upturned sides
- Removable for cleaning
- Compatible with DTMR & DTSR ranges only



Item	Description	Weight
SR	Reversible Shelves	5-10kg

Guard Rail

Features

- Stainless steel three-sided guard rail
- Increases the overall height of the trolley by approx. 60mm
- Compatible with DTSR range only
- Factory fitted



Item	Description	Weight
GR	Guard Rail	1.5kg

Clamping Infusion Stands - Trolley

Features

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trolleys, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum capacity on hooks 2kg

Overall Dimensions (w x d x h)

- 55 x 220 x 1200mm



Item	Description	Weight
DP/C	Clamping Infusion Stands - Trolley	1.8kg