



# General Purpose Trolleys

Stainless Steel Trolleys	186
Polymer Trolleys	188
Mild Steel Trolleys	189



## General Purpose Trolley Stainless Steel, Heavy Duty

### Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down (shelf size 920 x 535mm)
- Available with two or three shelves, shelf heights
  - Two shelf - 280 & 860mm
  - Three shelf - 220, 540 & 860mm
- Maximum capacity per shelf 75kg
- 125mm pressed steel castors

### Overall Dimensions (w x d x h)

- No Ring Buffers - 970 x 540 x 1020mm
- Ring Buffers - 1045 x 615 x 1020mm

### Standard Colours

- Stainless Steel

### Options

- C007 - 125mm Castors, Heavy Duty Pressed Steel, 4 x Non-Braking (set of 4)
- C008 - 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4)

Item	Description	Weight
WT/37/2	Stainless Steel - Heavy Duty, Two Shelves	28.0kg
WT/37/3	Stainless Steel - Heavy Duty, Three Shelves	37.5kg
WT/37/2/ZP010	Stainless Steel - Heavy Duty, Two Shelves, Ring Buffers	28.2kg
WT/37/3/ZP010	Stainless Steel - Heavy Duty, Three Shelves, Ring Buffers	37.7kg

## General Purpose Trolley Stainless Steel, Light Duty

### Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed liquid retaining sunken shelves (shelf size 540 x 430mm)
- Available with two or three shelves, shelf heights
  - Two shelf - 250 & 870mm
  - Three shelf - 250, 560 & 870mm
- Maximum capacity per shelf 50kg
- 100mm pressed steel castors

### Overall Dimensions (w x d x h)

- No Ring Buffers - 990 x 510 x 1035mm
- Ring Buffers - 1065 x 585 x 1035mm

### Standard Colours

- Stainless Steel



Item	Description	Weight
WT/33/SS/2	Stainless Steel - Light Duty, Two Shelves	16.0kg
WT/33/SS/3	Stainless Steel - Light Duty, Three Shelves	22.0kg
WT/33/SS/2/ZP005	Stainless Steel - Light Duty, Two Shelves, Ring Buffers	16.2kg
WT/33/SS/3/ZP005	Stainless Steel - Light Duty, Three Shelves, Ring Buffers	22.2kg

## General Purpose Trolley Polymer, Shallow Shelves

### Features

- Easy to clean corrosion resistant high density polyethylene shelves (not suitable for use with automated washers)
- 11mm raised edge retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria, mould, mildew & fungi
- Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two or three shelves in three sizes
  - Small - Shelf size 620 x 345mm
    - Easy grip handle & utility tray
    - Maximum capacity per shelf 60kg
    - 100mm pressed steel castors
  - Medium - Shelf size 690 x 475mm
    - Handle has built in recessed support channel to manage odd shaped items
    - Maximum capacity per shelf 60kg
    - 100mm pressed steel castors
  - Large - Shelf size 840 x 585mm
    - Handle has built in slot to hold scanners
    - Maximum capacity per shelf 75kg
    - 125mm pressed steel castors
- Self-assembly (supplied boxed)



### Standard Colours

- Uprights - Chrome
- Shelves - Blue

Item	Shelves	Weight
<b>Small - 800 x 465 x 900mm</b>		
8MY1627-24BU	2	12.9kg
8MY1627-34BU	3	15.4kg
<b>Medium - 873 x 595 x 900mm</b>		
8MY2030-24BU	2	14.9kg
8MY2030-34BU	3	18.8kg
<b>Large - 1022 x 703 x 936mm</b>		
8MY2636-25BU	2	18.1kg
8MY2636-35BU	3	23.1kg



8MY1627-34BU



8BC2030-3DMB

## General Purpose Trolley Polymer, Deep Shelves

### Features

- Easy to clean corrosion resistant high density polyethylene shelves (not suitable for use with automated washers)
- 70mm deep shelves retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria, mould, mildew & fungi
- Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two or three shelves in three sizes
  - Small - 800 x 546mm
  - Large - 948 x 686mm
- Max capacity shelf 68kg
- 100mm pressed steel castors
- Self-assembly (supplied boxed)

### Standard Colours

- Uprights - Chrome
- Shelves - Blue

Item	Shelves	Weight
<b>Small - 832 x 546 x 1041mm</b>		
8BC2030-2DMB	2	13.0kg
8BC2030-3DMB	3	17.0kg
<b>Large - 984 x 685 x 1041mm</b>		
8BC2636-2DMB	2	15.9kg
8BC2636-3DMB	3	21.1kg

## General Purpose Trolley Mild Steel, Laminate Shelves

### Features

- Laminate shelves
- Available with two or three shelves in two sizes
- Shelf heights
  - Two shelf - 210 & 845mm
  - Three shelf - 210, 525 & 845mm
- Shelf size
  - Small - 705 x 445mm
  - Large - 855 x 445mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel castors

### Standard Colours

- Paint - Grey
- Shelf - Spearmint

### Options

- ZP005 - Ring Buffer



GP015

Item	Description	Weight
<b>Small - 972 x 533 x 1008mm</b>		
GP005	Mild Steel - Two Laminate Shelves, Small	18.0kg
GP015	Mild Steel - Three Laminate Shelves, Small	23.0kg
<b>Large - 1122 x 533 x 1008mm</b>		
GP025	Mild Steel - Two Laminate Shelves, Large	20.0kg
GP035	Mild Steel - Three Laminate Shelves, Large	26.5kg

## General Purpose Trolley Mild Steel, Aluminium Trays

### Features

- Liquid retaining painted aluminium trays (depth 50mm)
- Available with two or three shelves in two sizes
- Shelf heights
  - Two shelf - 210 & 845mm
  - Three shelf - 210, 525 & 845mm
- Shelf size
  - Small - 705 x 445mm
  - Large - 855 x 445mm
- Maximum capacity per tray 40kg
- 100mm pressed steel castors

### Standard Colours

- Paint - Grey

### Options

- ZP005 - Ring Buffer



GP020

Item	Description	Weight
<b>Small - 972 x 533 x 1008mm</b>		
GP010	Mild Steel - Two Aluminium Trays, Small	16.5kg
GP020	Mild Steel - Three Aluminium Trays, Small	20.6kg
<b>Large - 1122 x 533 x 1008mm</b>		
GP030	Mild Steel - Two Aluminium Trays, Large	18.5kg
GP040	Mild Steel - Three Aluminium Trays, Large	23.7kg