



## Portering Chairs

Reer Steer, Sliding Foot Rest	264
Front Steer, Sliding Foot Rest	265
Front Steer, Hinged Foot Rests	266
Bariatric, Reer Steer, Sliding Foot Rest	267
Options	266



CE



G/200/RS/VBB

### Portering Chair Rear Steer, Nesting, Sliding Foot Rest

#### Features

- Nests to occupy reduced floor space and encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should patient stand on the step
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 200mm fixed wheels/150mm total stop braking castors (cushioned)
- Maximum capacity 190.0kg

#### Dimensions (w x d x h)

- Overall Dimensions - 660 x 840 x 1040mm
- Seat - 540 x 420mm at Height 530mm
- Back - 490 x 455mm

#### Standard Colours

- Paint - Grey
- Upholstery - Bristol Blue

#### Options

- CYHP/D1 - Cylinder Holder, CD/D
- CYHP/E1 - Cylinder Holder, E
- G/CLS - Calf/Leg Supports (pair)
- G/DP - Infusion Pole
- G/SB - Lap Seatbelt, Factory Fitted
- G/CARD/LOCK - Card/Coin Lock, Factory Fitted
- G/CARD/STATION - Head of Station
- 2LOKE011 - Release Card



X-Ray/Notes Pocket

Item	Description	Weight
G/200/RS/VBB	Rear Steer, Nesting, Sliding Foot Rest	29.0kg

Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer, Nesting, Sliding Foot Rest in Bristol Blue - G/200/RS/VBB



VBB VAQ VFR VGH VSD VOB VLI VCH VCB VPE



CE

### Portering Chair Front Steer, Nesting, Sliding Foot Rest

#### Features

- Nests to occupy reduced floor space and encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should patient stand on the step
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg

#### Dimensions (w x d x h)

- Overall Dimensions - 660 x 840 x 1040mm
- Seat - 540 x 420mm at Height 530mm
- Back - 490 x 455mm

#### Standard Colours

- Paint - Grey
- Upholstery - Bristol Blue

#### Options

- CYHP/D1 - Cylinder Holder, CD/D
- CYHP/E1 - Cylinder Holder, E
- G/CLS - Calf/Leg Supports (pair)
- G/DP - Infusion Pole
- G/SB - Lap Seatbelt, Factory Fitted
- G/CARD/LOCK - Card/Coin Lock, Factory Fitted
- G/CARD/STATION - Head of Station
- 2LOKE011 - Release Card



X-Ray/Notes Pocket

Item	Description	Weight
G/305/FS/VBB	Front Steer, Nesting, Sliding Foot Rest	29.0kg

Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Sliding Foot Rest in Lime - G/305/FS/VLI



VBB VAQ VFR VGH VSD VOB VLI VCH VCB VPE

G/305/FS/VBB





CE



G/205/FS/VBB

### Portering Chair Front Steer, Nesting, Hinged Foot Rests

#### Features

- Nests to occupy reduced floor space and encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Folding/hinging foot rests aids access as well as supporting feet during transit
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg

#### Dimensions (w x d x h)

- Overall Dimensions - 685 x 710 x 1040mm
- Seat - 540 x 420mm at Height 530mm
- Back - 490 x 455mm

#### Standard Colours

- Paint - Grey
- Upholstery - Bristol Blue

#### Options

- CYHP/D1 - Cylinder Holder, CD/D
- CYHP/E1 - Cylinder Holder, E
- G/CLS - Calf/Leg Supports (pair)
- G/DP - Infusion Pole
- G/SB - Lap Seatbelt, Factory Fitted
- G/CARD/LOCK - Card/Coin Lock, Factory Fitted
- G/CARD/STATION - Head of Station
- 2LOKE011 - Release Card



X-Ray/Notes Pocket

Item	Description	Weight
G/205/FS/VBB	Front Steer, Nesting, Hinged Foot Rests	27.5kg

Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Hinged Foot Rests in Flame Red - G/205/FS/VFR



### Portering Chair Bariatric, Rear Steer, Sliding Foot Rest



CE



G/500/RS/VBB

#### Features

- Heavy duty bariatric portering chair
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should patient stand on the step
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles mounted to rear uprights
- Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg
- Maximum capacity 380.0kg

#### Dimensions (w x d x h)

- Overall Dimensions - 930 x 1000 x 1000mm
- Seat - 740 x 670mm at Height 530mm
- Back - 700 x 480mm

#### Standard Colours

- Paint - Grey
- Upholstery - Bristol Blue

#### Options

- G/500/RS/CYL - Cylinder Holder, CD/D
- G/CLS/500 - Calf/Leg Supports (pair)
- G/DP - Infusion Pole
- G/CARD/LOCK - Card/Coin Lock, Factory Fitted
- G/CARD/STATION - Head of Station
- 2LOKE011 - Release Card



X-Ray/Notes Pocket

Item	Description	Weight
G/500/RS/VBB	Bariatric, Rear Steer, Sliding Foot Rest	43.8kg

Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Bariatric, Rear Steer, Sliding Foot Rest in Coastal Blue - G/500/RS/VCB



## Cylinder Holders

### Features

- Compatible with 100mm diameter cylinders (C/CD/D/E)
- Mounts to IV support tube at rear of chair using a handwheel
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models

### Standard Colours

- Paint - Grey



Item	Description	Weight
CYHP/D1	C/CD/D Size	1.0kg
CYHP/E1	E Size	1.8kg
G/500/RS/CYL	Bariatric - C/CD/D Size	1.0kg

## Calf / Leg Support - Pair

### Features

- Moulded calf/leg support on captive retractable support tube
- Retracts beneath the seat when not in use
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Supplied as a pair

### Standard Colours

- Stainless Steel/Black



Item	Description	Weight
G/CLS	Calf/Leg Support, Pair	2.3kg
G/CLS/500	Bariatric Calf/Leg Support, Pair	2.3kg

## Infusion Pole

### Features

- Stainless steel infusion pole with a U shaped hook (length of pole 980mm)
- Mounts to IV support tube at rear of chair using a handwheel
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models



Item	Description	Weight
G/DP	IV Pole	1.1kg

## Lap Seatbelt

### Features

- Lap seat belt with snap buckle
- Compatible with G/200/RS, G/205/FS & G/305/FS models
- Factory fitted

### Standard Colours

- Black



Item	Description	Weight
G/SB	Lap Seatbelt	0.4kg

## Card/Coin Lock

### Features

- Compatible with old & new one pound coin (euro compatible)
- Coin/card remains in the lock until the chair is connected to another chair or the head of station
- Plastic release card available for use by porters
- Factory fitted
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Replacement card/coin lock cassette available (requires cassette removal tool)



Item	Description	Weight
G/CARD/LOCK	Card/Coin Lock, Factory Fitted	1.0kg
2LOKE011	Plastic Release Card	0.01kg
2LOCK217	Card/Coin Lock Cassette, New £1 Coin Compatible	0.2kg
2LOCK218	Cassette Removal Tool	0.01kg

## Card/Coin Lock Wall Bracket

### Features

- Head of station for chairs fitted with card/coin locks
- Requires fitting to the fabric of the building (wall fixings not supplied)
- Chain length 1000mm
- One required for each row of chairs, (G/200/RS, G/205/FS & G/305/FS)
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models



Item	Description	Weight
G/CARD/STATION	Head of Station	1.6kg