



Procedure Trolleys & Carts

Plaster Trolleys 272

Mobile Workstation 273

Procedure Carts 274



WT/96/3/SS/4BK
(fitted with optional flap,
louver panel & ring buffers)

2PLCO024
optional lin bin

WT/96/2/SS/2BK

WT/96/1/SS

Plaster Trolleys - Stainless Steel

Features

- Seamlessly welded easy clean stainless steel plaster trolley (304 grade)
- Available with one/two or three drawers
- Integral bucket ring (supplied with a 15 litre stainless steel bucket)
- 100mm pressed steel anti-static non marking castors

Dimensions (w x d x h)

- Overall Dimensions - 1150 x 540 x 915mm
- Drawers - 650 x 410 x 100mm (internal)
- Work Surface - 745 x 515 at height 915mm
- Lower Shelf - 746 x 515 at height 175mm

Options

- WT/96/FLAP - Lift Up Flap
- WT/96/SS/LP - Louver Panel
- 2PLCO024 - Plastic Lin Bin
- ZP010 - Ring Buffers
- 2BK - 2 x Braking Castors / 2 x Non-Braking Castors
- 4BK - 4 x Braking Castors

Standard Colours

- Stainless

Item	Description	Weight
WT/96/1/SS	One Drawer & Bucket	33.0kg
WT/96/2/SS	Two Drawers & Bucket	42.0kg
WT/96/3/SS	Three Drawers & Bucket	51.0kg

Mobile Workstation - Mild Steel, Low Level

Features

- Low level mobile workstation constructed from mild steel
- Recessed liquid retaining top
- Single point locking
- To maximise safety, the unit is fitted with an anti-tilt mechanism, allowing only one drawer to be opened at any one time
- Full extension drawers, heights 85, 170 & 255mm (max load 30.0kg)
- Two tone colour keyway, drawer fronts available in a range of colours
- Low level corner bumpers to protect the fabric of the building
- 80mm heavy duty castors, 2 off braking

Dimensions (w x d x h)

- Overall - 615 x 500 x 810mm
- Drawers - 480 x 375 x 85/170/255mm (internal)

Standard Colours

- Body - Light Grey
- Drawers - Dark Grey

Item	Description	Weight
MW005	2 x 85mm, 1 x 170mm, 1 x 255mm Drawers, Cam Lock	49.0kg





Procedure Carts - Polymer, Metro Flexline™

Features

- Range of pre-configured carts to suit a wide range of applications
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Upper work surface provides a writing surface/space for preparation work
- Three locking options (CAM/electronic locking & security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152, 229mm (max load 24kg)
- Incorporates Microban offering antimicrobial protection
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm castors, 2 x Directional Castors / 2 x Braking Castors

Dimensions (w x d x h)

- Drawers - 514 x 432 x 76/152/229mm (internal)

Standard Colours

- Body - Light Taupe
- Drawers Handles - Model Dependant

Options

- 8FL190 - Drawer Label Holders (x 10)



Drawer Label Holders

This is a small selection of a large configurable range. Please contact our customer services team to discuss your requirements

Item	Description	Overall Dimensions (w x d x h)	Weight
8FLCRP4	Emergency Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Code Blue Handles, Security Seal Lock	1056 x 607 x 1273mm	91.0kg
8FLPED1	Paediatric Emergency Cart - 8 x 76mm, 1 x 152mm Drawers, Coloured Handles, Security Seal Lock	1056 x 607 x 1273mm	91.0kg
8FLTMNT2	Treatment Cart - 2 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles, CAM Lock	819 x 568 x 1064mm	51.0kg
8FLCOM-LT	Laptop Cart - 1 x Keyboard Tray, 1 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles, Keyless Entry Touch Pad	819 x 568 x 1064mm	83.0kg
8FLCLINE	Central Line Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, Keyless Entry Touch Pad	819 x 568 x 1064mm	91.0kg
8FLCCU1	Critical Care Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Red Handles, Keyless Entry Touch Pad	840 x 569 x 1807mm	90.0kg



Fitted To Trolley

	Emergency 8FLCRP4	Paediatric Emergency 8FLPED1	Treatment 8FLTMNT2	Laptop 8FLCOM-LT	Central Line 8FLCLINE	Critical Care 8FLCCU1
76mm Drawer Divider Kit			2		2	
152/229mm Drawer Divider Kit			2		1	
76mm Drawer Divider Tray & Dividers	3					
152mm Drawer Divider Tray & Dividers	2					
Side Bin, Locking	4				2	4
Side Bin, Non-Locking				3	2	4
Waste Basket (28 Litre)					✓	
Oxygen Cylinder Holder	✓	✓				
Lockable Sharps Container					✓	
Cord Manager	✓					
Adjustable Defibrillator Tray	✓	✓				
Hospital Grade 6, Outlet Strip & Holder	✓					
Suction Pump Shelf	✓	✓				
Rear Fitting Backboard	✓	✓				
Laptop Arm				✓		
Pull Out Side Shelf	✓	✓				
Two Hook Infusion Pole	✓	✓				
Plastic Security Seals, (x 100)	✓	✓				
Drawer Seal Lock Bar		✓				
Lock Alert VI – Auto Lock				✓	✓	✓
Overbridge with Top Shelf & Hanger Rail						✓
Six Compartment Tilt Bin for Overbridge						✓



Procedure Carts - Polymer, Metro Starsys™

Features

- Range of pre-configured carts to suit a wide range of applications
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Upper work surface provides a writing surface/space for preparation work
- Five locking options, (electronic touchpad/pushbutton/keyless, CAM & security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152, 229mm (max load 11/18kg)
- Incorporates Microban offering antimicrobial protection
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm castors, 2 x Directional Castors / 2 x Braking Castors

Standard Colours

- Body - Light Taupe
- Drawers Handles - Model Dependant

Options

- 8SXLABKIT - Drawer/Door Label Holders, (x 10)



Drawer / Door Label Holders

This is a small selection of a large configurable range. Please contact our customer services team to discuss your requirements

Item	Description	Overall Dimensions (w x d x h)	Weight
8SXRANES	Anaesthesia Cart - 4 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock	882 x 633 x 1136mm	104.0kg
8SXRSDRS	Dressing Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock	577 x 633 x 1060mm	83.0kg
8SXRSTRAMA	Trauma Cart - 3 x 76mm, 4 x 152mm Drawers, Slate Blue Handles, CAM Lock	882 x 633 x 1136mm	116.0kg
8SXRSCCU	Critical Care Cart - 1 x 76mm, 1 x 152mm, 3 x 229mm Drawers, Slate Blue Handles	882 x 633 x 1212mm	107.0kg
8SXRARTH	Arthroscopic Cart - 6 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock	882 x 633 x 1136mm	110.0kg
8SXRIV	IV Therapy Cart - 5 x 76mm, 2 x 229mm Drawers, Green Handles, CAM Lock	882 x 633 x 1136mm	106.0kg
8SXRSCAST	Cast Cart - 2 x 76mm, 2 x 152mm, 2 x 229mm Drawers, Slate Blue Handles	882 x 633 x 1212mm	109.0kg



Fitted To Trolley

	Anaesthesia 8SXRANES	Trauma 8SXRSTRAMA	Critical Care 8SXRSCCU	Arthroscopic 8SXRARTH	IV Therapy 8SXRIV	Cast 8SXRSCAST
76mm Drawer Divider Kit		2	1		4	1
Left Swing Out Pod, Three Tilt Bins			✓			✓
Left Swing Out Pod, Three Locking Tilt Bins (CAM Lock)					✓	
Left Swing Out Pod, Top Locking Tilt Bin (CAM Lock)						
Pod Tray With Dividers		✓		✓		
Waste Basket & Holder		✓	✓	✓		✓
Right Fixed Pod, One Tilt Bin			✓			✓
Right Swing Out Pod				✓		
Right Swing Out Pod, Three Locking Tilt Bins (CAM Lock)					✓	
Right Swing Out Pod, Top Locking Tilt Bin (CAM Lock)		✓				
Tilt Bin				✓		
Middle Non Locking Tilt Bin		✓				
Sharps Container Bracket		✓		✓		
Left Swing Out Locking Pod, Three Tilt Bins (CAM Lock)	✓					
Right Swing Out Locking Pod, Three Tilt Bins (CAM Lock)	✓					
Overbridge with One Shelf & Accessory Rail		✓	✓			
Overbridge with Two Accessory Rails	✓			✓		✓
Two Hook Infusion Pole		✓			✓	
Wire Basket		✓	✓	✓		✓
Universal Clamp	✓					
Laryngoscope Blade/Specimen Holder	✓					
Half Size Utility Bin (Metal)	✓					
Full Size Utility Bin (Metal)	✓					
Label/Tape Dispenser (x 3)	✓					
Plastic Utility Bins with Clear Covers (x 6)	✓		✓	✓		



Procedure Carts - Polymer, Metro Basix™

Features

- Lightweight metal procedure carts
- Recessed liquid retaining plastic work surface
- Three locking options
 - CAM (Compact)
 - Push button (Trauma)
 - Security seal compatible (Anaesthesia)
- Full extension drawers, depths 76, 152 & 230mm
 - Compact - Max load 7kg
 - Trauma/Anaesthesia - Max load 24kg
- Drawer fronts available in a range of colours
- 125mm castors, total lock at front & directional lock at rear

Dimensions (w x d x h)

- Compact - 489 x 470 x 876mm
- Drawers (internal)
 - 76mm - 381 x 432 x 76mm
 - 152mm - 381 x 432 x 152mm
 - 230mm - 381 x 432 x 230mm
- Trauma/Anaesthesia Cart - 762 x 613 x 1440mm
- Drawers (internal)
 - 76mm - 432 x 565 x 76mm
 - 152mm - 432 x 565 x 152mm
 - 230mm - 432 x 565 x 230mm

Standard Colours

- Body - Light Taupe
- Drawers - Dark Blue
(see below for alternative colours)




Label Holders

Options

- 8MBA172 - Label Holders (x 12)

Alternative drawer colours are available. Change suffix 'DB' on the product code with any of the following i.e. - Basix Cart with red drawers 8MBC1210RE



Item	Description	Weight
8MBC1210DB	Compact - 1 x 76mm, 2 x 152mm, 1 x 230mm Drawers	63.0kg
8MBX3210TRAUMA	Trauma - 3 x 76mm, 2 x 152mm, 1 x 230mm Drawers	78.0kg
8MBX1310ANES4	Anaesthesia - 1 x 76mm, 3 x 152mm, 1 x 230mm Drawers	81.0kg



8MBC1210YE



8MBX3210TRAUMARE



8MBX1310ANES4TE