



## Resuscitation Trolleys & Carts

Mild Steel Resuscitation Trolleys 282

Polymer Resuscitation Carts 284

Airways / Difficult Intubation Trolley 287

Compact Resuscitation Trolleys 288



WT210

WT200

WT250/CYL/D/E  
shown with optional cylinder holder WT/CYL/D/E

## Resuscitation Trolleys - Mild Steel

### Features

- Popular resuscitation trolley
- Drawer contents are secured when the lower cupboard is closed (security seal compatible)
- Supplied with a plastic drawer liner with dividers
- Recessed liquid retaining stainless steel work surface (usable 630 x 445)
- Lift up work flap (usable 305 x 455mm)
- Two hook infusion pole mounted to rear of trolley
- Low level bumper bar
- Push handles one end
- 125mm pressed steel castors

### Dimensions (w x d x h)

- Drawers (internal)
  - 70mm - 600 x 430 x 70mm
  - 96mm - 600 x 430 x 96mm
  - 155mm - 600 x 430 x 155mm
- Lower Cupboard (internal)
  - 625 x 475 x 295mm

### Standard Colours

- Paint - Red
- Lift Up Flap - White

Drop down cupboard door prevents drawers from being opened. Can be secured using security seals.

| Item  | Description                                       | Overall Dimensions (w x d x h) | Weight |
|-------|---|--------------------------------|--------|
| WT200 | 2 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 999mm              | 74.0kg |
| WT210 | 3 x 96mm Deep Drawers, Lower Cupboard             | 825 x 575 x 999mm              | 75.1kg |
| WT250 | 4 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 1165mm             | 86.0kg |

## Cardiac Board

### Features

- Laminated cardiac board approx 500 x 400 x 17mm, complete with fixing brackets
- Positioned at the rear of the trolley for easy access
- Factory fitted

| Item      | Description   | Weight |
|-----------|---------------|--------|
| WT/CBOARD | Cardiac Board | 1.5kg  |



## Cylinder Holders

### Features

- WT/CYL/D/E compatible with AD, CD, D, DD, PD, RD, ZD & E cylinders, max diameter 100mm
- WT/CYL/F compatible with F, LF, VF, HX & ZX cylinders, max diameter 145mm
- Maximum two cylinder holders per trolley
- Factory fitted

| Item       | Description                  | Weight |
|------------|------------------------------|--------|
| WT/CYL/D/E | D/E Size, Max Diameter 100mm | 3.0kg  |
| WT/CYL/F   | F Size, Max Diameter 145mm   | 3.4kg  |



## Drawer Tidy

### Features

- Plastic drawer liner with removable dividers
- One tidy supplied as standard with each trolley

| Item  | Description | Weight |
|-------|-------------|--------|
| WT/DT | Drawer Tidy | 1.6kg  |



## Four Gang Socket

### Features

- Switched & fused four gang socket
- Supplied with 3 metres of beige cable & fitted with a UK 3 pin plug
- Fitted to the rear of the trolley, complete with a cleat to secure the cable when not in use
- Factory fitted



| Item       | Description      | Weight |
|------------|------------------|--------|
| E/SOCKET/4 | Four Gang Socket | 1.4kg  |

## Guard Rail

### Features

- Stainless steel guard rail
- Increases the height of the trolley by 60mm
- Compatible with WT200, WT210 and WT250 models
- Factory fitted



| Item     | Description | Weight |
|----------|-------------|--------|
| G/RAIL-C | Guard Rail  | 1.5kg  |

## Security Seals

### Features

- Snap on plastic security
- Pack of 100



| Item     | Description    | Weight |
|----------|----------------|--------|
| 2PLCO026 | Security Seals |        |

## Heavy Duty Castors

### Features

- 125mm heavy duty pressed steel castors
- Offers additional cushioning, lowering noise levels
- Supplied as a set of 4 castors
  - 2 x Braking, 2 x Non-Braking
  - 4 x Non-Braking
- Factory fitted

| Item | Description                             | Weight |
|------|---|--------|
| C011 | 4 x Non-Braking (set of 4)              | 4.0kg  |
| C012 | 2 x Braking, 2 x Non-Braking (set of 4) | 5.0kg  |



8LECCRP4  
(defibrillator not included)



8LECCRP2  
(defibrillator not included)



Drawer Label Holders

**Features**

- Range of pre-configured emergency carts
- Incorporates Microban offering antimicrobial protection
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Single seal mechanism secures the top compartment, drawers and tilt out side bins (security seal compatible)
- Top storage compartment with clear removable cover provides instant access to first line equipment
- Strapless defibrillator platform on swing arm adjusts to accommodate various models
- Full extension lightweight polymer drawers, depths 76, 152, 229 & 305mm, (max load 24kg)
- Ergonomic handle with fifth wheel steering provides maximum control when attending an emergency situation
- Two tone colour keyway
- 127mm castors, 2 x Directional / 2 x Braking

**Dimensions (w x d x h)**

- Drawers - 514 x 432 x 76 / 152 / 229 / 305mm (internal)

**Standard Colours**

- Body - Code Blue
- Drawers - Light Taupe/Code Blue

**Options**

- 8FL190 - Drawer Label Holders (x 10)

| Item     | Description                            | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| 8LECCRP2 | 1 x 76mm, 3 x 152mm Drawers            | 965 x 575 x 984mm              | 84.0kg |
| 8LECCRP3 | 2 x 76mm, 1 x 152mm, 1 x 229mm Drawers | 965 x 575 x 984mm              | 88.0kg |
| 8LECCRP4 | 2 x 76mm, 1 x 152mm, 1 x 305mm Drawers | 965 x 575 x 1070mm             | 87.0kg |



defibrillator not included

|                                |   |   |   |
|--------------------------------|---|---|---|
| Top Cavity Tray                |   | ✓ |   |
| Locking Side Bin               |   | 3 |   |
| Waste Basket (28 Litre)        | ✓ |   |   |
| Lockable Sharps Container      | ✓ |   |   |
| Adjustable Defibrillator Tray  | ✓ |   | ✓ |
| Suction Pump Shelf             |   | ✓ |   |
| Backboard (Front Fitting)      |   | ✓ | ✓ |
| Backboard (Rear Fitting)       | ✓ |   |   |
| Oxygen Cylinder Holder         | ✓ | ✓ | ✓ |
| Two Hook Infusion Pole         | ✓ | ✓ | ✓ |
| Storage Bin                    | ✓ | ✓ | ✓ |
| Plastic Security Seals (x 100) | ✓ | ✓ | ✓ |



8FL21P3369  
(fitted with optional two hook infusion pole, oxygen cylinder holder & defibrillator platform. Defibrillator not included)

## Resuscitation Carts - Polymer, Metro Flexline™

### Features

- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban® offering antimicrobial protection
- Upper work surface provides a writing surface/space for preparation work
- Single seal mechanism secures the drawers (security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152, 229mm (max load 24kg)
- Two tone colour keyway
- 127mm castors, 2 x Directional / 2 x Braking

### Dimensions (w x d x h)

- Overall Dimensions - 819 x 568 x 914mm
- Drawers - 514 x 432 x 76 / 152 / 229mm (internal)

### Standard Colours

- Body - Light Taupe
- Drawers - Light Taupe/Red

### Options

- 8FL141 - 76mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL142 - 152mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL315 - Two Hook Infusion Pole
- 8FL234 - 100mm Diameter Oxygen Cylinder Holder
- 8FL190 - Drawer Label Holders (x 10)
- 8LEC320 - Plastic Security Seals (x 100)
- 8LEC304 - Defibrillator Platform

| Item       | Description  | Weight |
|------------|--|--------|
| 8FL21P3369 | 2 x 76mm, 1 x 152mm, 1 x 229mm Drawers, Red Drawer Pulls, Security Seal Lock | 68.0kg |

## Airways / Difficult Intubation Trolley - Mild Steel

### Features

- Designed specifically for airways / difficult intubation in conjunction with resuscitation officers
- Recessed liquid retaining top
- Available as non-locking or fitted with a side mounted tamper evident panel to secure the contents of the drawers (security seal compatible)
- 135mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Two accessory rails fitted to L/H side
- Low level corner bumpers to protect the fabric of the building
- 80mm heavy duty castors, 2 off braking

### Dimensions (w x d x h)

- Overall - 615 x 500 x 915mm
- Drawers - 480 x 375 x 135mm (internal)

### Standard Colours

- Body - White
- Drawers - White

### Options



CT/IV

CT/BIN

SBH/15

CT/AGH

GMBH/15

CT/DH/A4

CT/CYD

CT/AIR

CT/CR5 & CT/CR2



MW015 Tamper Evident Panel



MW010

| Item  | Description                             | Weight |
|-------|---|--------|
| MW010 | 5 x 125mm Drawers, Non Locking          | 55.0kg |
| MW015 | 5 x 125mm Drawers, Tamper Evident Panel | 57.0kg |



RTM1/10/TD/3S2D

## Resuscitation Trolleys - Compact

### Features

- Single column standard level resuscitation trolleys
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trolley to support a wide range of accessories (see options)
- Fitted with a back panel & tamper evident door providing security for the contents of the trays (security seal compatible)
- Trolley pre-configured with the following...
  - Cylinder Holder - Compatible with C/CD/D size cylinders
  - Single Glove / Mask Box Holder - Accommodates glove & mask boxes up to 140 mm x 95 mm
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with six sheets of perforated coloured label strips - white, blue, green, yellow & red (not shown on trolley)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm polymer castors offering smooth surfaces to assist cleaning, 2 off braking

### Dimensions (w x d x h)

- Overall Dimensions - 550 x 550 x 993mm
- Overall Dimensions (with overbridge) - 600 x 550 x 1620mm
- Plastic Top
  - Single Column - 455 x 455 (internal) at height 968mm
- Trays
  - Shallow - 355 x 405 x 100mm (internal)
  - Deep - 355 x 405 x 150mm (internal)

### Standard Colours

- Paint - Red

| Item            | Description                | Weight |
|-----------------|----------------------------|--------|
| RTM1/10/TD/3S2D | 3 x 100mm, 2 x 150mm Trays | 33.3kg |
| RTM1/10/TD/4S1D | 4 x 100mm, 1 x 150mm Trays | 32.9kg |
| RTM1/10/TD/6S   | 6 x 100mm Trays            | 34.0kg |

## Overbridge

### Features

- Fixes onto the rear uprights of the trolley
- Provides a fixing location for tilt bins (order separately)
- Factory fitted

### Standard Colours

- Paint - Grey White

RTM1/10/TD/3S2D  
shown with Overbridge & Tilt Bins

| Item     | Description            | Weight |
|----------|------------------------|--------|
| CT/OB/MS | Overbridge, Mild Steel | 2.0kg  |

## Tilt Bins

### Features

- Available with 4, 5 or 6 tilting compartments
- Clear tilt bins provide visible storage & assists organisation
- Compartments tilt forward to a 45° angle allowing easy access of contents & can be removed to aid refilling/cleaning
- Label holder for easy identification of the contents
- Retaining bar secures compartments in place, preventing bins from tipping forward during transit
- Attaches to an overbridge (order separately)

### Dimensions (w x d x h)

- Compartment Dimensions
- TBIN4/MS - 143 x 127 x 190mm
  - TBIN5/MS - 108 x 95 x 146mm
  - TBIN6/MS - 89 x 67 x 95mm

### Standard Colours

- Plastic - Cream
- Paint - Grey White



45° Angle Tilt



Label Holder

Top: TBIN6/MS  
Middle: TBIN5/MS  
Bottom: TBIN4/MS

| Item     | Description      | Dimensions (w x d x h) | Weight |
|----------|------------------|------------------------|--------|
| TBIN4/MS | Four Compartment | 600 x 171 x 207mm      | 2.3kg  |
| TBIN5/MS | Five Compartment | 600 x 135 x 164mm      | 1.6kg  |
| TBIN6/MS | Six Compartment  | 600 x 93 x 113mm       | 0.9kg  |

### Options (refer to page 66 for additional information)

- CT/AGH - Alcohol Gel Holder
- CT/AIR - Airway Tube Holder
- CT/BIN - Debris Bin & Bracket
- CT/DH/A4 - Document Holder, A4
- CT/CR2 - Catheter Rack, Two Section
- CT/CR5 - Catheter Rack, Five Section
- CT/FS - Flat Shelf - Removable
- CT/IV - Infusion Pole & Bracket
- CT/OB/MS - Overbridge
- TBIN4/MS - Tilt Bin, 4 Compartment (overbridge required)
- TBIN5/MS - Tilt Bin, 5 Compartment (overbridge required)
- TBIN6/MS - Tilt Bin, 6 Compartment (overbridge required)
- CT/WF - Work Flap
- GMBH/15 - Glove / Mask Box Holder - One Box, 15mm Hook
- SBH/15 - Sharps Box Holder - One Box, 15mm Hook
- CT/SBB/15 - Sharps Box Bracket - Universal (Daniels / Frontier / Sharpak)
- 2PLCO026 - Security Seals, Pack of 100