

METTI
(Muallem ET
Tube Introducer)

for difficult intubation
and tube exchange

- core made of plastic,
semi-rigid (flexible)
- flexible, pre-formed tip
- graduation marks
- for single use, sterile



**INTUBATION-
CATHETER**

for fiberoptically assisted
intubation ($\leq 4,0$ mm)

- lumen-design,
O₂-Connector included
- graduation marks
- for single use, sterile



TUBE EXCHANGER

for extubation and
exchange of Endo-
tracheal Tubes

- lumen-design,
O₂-Connector included
- graduation marks
- for single use, sterile



STYLETS

malleable intubation
stylet

- soft, atraumatic tip
- with and without
moveable
silicone connector
- autoclavable



ORDER INFORMATION

METTI

Size	Length	for ET Tube	REF	Box
12 CH	80 cm	≥ 5.0 mm	33-12-800-1	5
14 CH	80 cm	≥ 6.0 mm	33-14-800-1	5

Intubation-Catheter

Size	Length	for ET Tube	for Fiberoptic	REF	Box
19 CH	56 cm	≥ 7.0 mm	≤ 4.0 mm	33-60-560-1	5

Tube Exchanger

Size	Length	for ET Tube	REF	Box
11 CH	80 cm	≥ 4.0 mm	33-45-800-1	5
14 CH	80 cm	≥ 5.0 mm	33-46-800-1	5
19 CH	80 cm	≥ 7.0 mm	33-47-800-1	5

Stylet

Size	Length	for ET Tube	without silicone connector	Box	with moveable silicone connector (15 mm I.D.)	
					REF	Box
small	34 cm	2.5 - 3.5 mm	REF 90-10-111	25	REF 90-10-001	10
medium	38 cm	4.0 - 4.5 mm	REF 90-10-222	25	REF 90-10-002	10
large	45 cm	5.0 - 6.5 mm	REF 90-10-333	25	REF 90-10-003	10
large	45 cm	7.0 - 11.0 mm	REF 90-10-444	25	REF 90-10-004	10